



## New Arrival | Cloud-Clone Recombinant Antibody


With the development of genetic engineering, the technology of monoclonal antibody preparation got rid of the limitation of traditional hybridoma technology, which resolved a series of ethics problems caused by tumor mice. Compared with traditional preparation technology, in aspects of species choice, antibody screening, titer control and mass production, genetic engineering antibodies show pretty advantages. Cloud-Clone Corp. provides one stop preparation service of full-length genetic engineering monoclonal antibody for multiple species such as mouse, rat, rabbit, etc.




### Flow Cytometry



Human peripheral blood lymphocytes were stained with purified CD4 (red) or purified mouse IgG1,  $\kappa$  isotype control (green), followed by anti-mouse IgG FITC.




Human peripheral blood lymphocytes were stained with CD8-APC (green) or mouse IgG1,  $\kappa$  isotype control-APC (blue).




Human peripheral blood lymphocytes were stained with purified CD45 (red) or purified mouse IgG1,  $\kappa$  isotype control (green), followed by anti-mouse IgG FITC.


## Immunohistochemical detection



IGF1: DAB staining on IHC-P;  
Samples: Human Kidney Tissue;




SAA: DAB staining on IHC-P;  
Samples: Human Kidney Tissue;




Cys-C: DAB staining on IHC-P;  
Samples: Human Kidney Tissue;


## Construction of immunoassay kits




Specificity and Linear Range of  
PINP Recombinant Antibody Pairs



Comparison of PINP Recombinant  
Antibody Pairs with Similar Products



Specificity and Linear Range of  
CRP Recombinant Antibody Pairs



Comparison of CRP Recombinant  
Antibody Pairs with Similar Products

## Cloud-Clone Recombinant Antibody:

Product No.	Host	Product Name	Product No.	Host	Product Name
RAB167Hu21	Mouse	Recombinant Antibody to Cluster Of Differentiation 4 (CD4)	RAA884Hu21	Mouse	Recombinant Antibody to Cluster Of Differentiation 26 (CD26)
RAA519Hu21	Mouse	Recombinant Antibody to Apolipoprotein A1 (APOA1)	RAA050Hu21	Mouse	Recombinant Antibody to Insulin Like Growth Factor 1 (IGF1)
RAA087Ca21	Mouse	Recombinant Antibody to Monocyte Chemotactic Protein 1 (MCP1)	RAA478Hu21	Mouse	Recombinant Antibody to Cardiac Troponin I (cTnI)
RAB167Mu31	Mouse	Recombinant Antibody to Cluster Of Differentiation 4 (CD4)	RAA896Hu21	Mouse	Recombinant Antibody to Cystatin C (Cys-C)
RAB167Mu32	Rat	Recombinant Antibody to Cluster Of Differentiation 4 (CD4)	RAA957Hu21	Mouse	Recombinant Antibody to Procollagen I N-Terminal Propeptide (PINP)
RAA821Hu21	Mouse	Recombinant Antibody to C Reactive Protein (CRP)	RAB030Hu21	Mouse	Recombinant Antibody to Protein Tyrosine Phosphatase Receptor Type C (CD45)
RAA821Hu22	Mouse	Recombinant Antibody to C Reactive Protein (CRP)	RAQ127Hu21	Mouse	Recombinant Antibody to Cluster Of Differentiation 8b (CD8b)
RAA821Hu23	Mouse	Recombinant Antibody to C Reactive Protein (CRP)	RAB837Hu21	Human	Recombinant Antibody to Cluster Of Differentiation 25 (CD25)
RAA885Hu21	Mouse	Recombinant Antibody to Serum Amyloid A (SAA)	RAD117Hu21	Mouse	Recombinant Antibody to Cluster Of Differentiation 3e (CD3e)

## Cloud-Clone Corp.

Original manufacturer of assay reagent, analysis equipment and animal experiment service.

23603 W. Fernhurst Dr., Unit 2201, Katy, TX 77494, USA Tel: 001-832-538-0970 Toll free: 888-960-7402 (In the USA)

[www.cloud-clone.com](http://www.cloud-clone.com) [www.cloud-clone.us](http://www.cloud-clone.us)